



Vessel Cannula

ORDER CODE: CM-4303

Vein graft vessel cannulas are used with heart, cardiopulmonary bypass, or vascular surgery to help check for leaks in a harvested vein that will be used in a graft. The cannulas feature an inlet female threaded Luer connector. The outlet tip is blunted for patient safety.

Outlet Tip Outer Diameter: Approximately 2mm

The construction is non-DEPH and of class 6 Polyvinyl material.

Pouch Contents:

Vessel Cannula
Luer Cap

EO Sterilized & Packaged 25/Box



Quality Management System

Registered to ISO 13485:2003 under CMDCAS
Annex V of MDD 93/42/EEC, FDA

CardioMed Supplies Inc.
199 Saint David Street
Lindsay, Ontario, Canada
K9V 5K7

1-800-387-9757
www.CardioMed.com