

Embolectomy Catheter

Maximum ease, minimum trauma

CardioMed's Embolectomy catheters facilitate rapid removal of emboli and thrombi from the arterial system.

Reliable and consistent balloon symmetry provides uniform contact with vessel walls

Rounded, soft tip enables easy insertion for minimized trauma

Catheter length and balloon inflation size are conveniently printed on proximal end for rapid identification

All embolectomy catheters come with a stylet and 3cc inflation syringe



Embolectomy Catheter

Ordering Information

Item	Description	French Size	Catheter Length cm	Maximum Inflation (Liquid) cc	Balloon Diameter mm (+/- 2mm)
CM702-40	2F x 40 cm Embolectomy Catheter	2	40	0.15	4
CM702-60	2F x 60 cm Embolectomy Catheter	2	60	0.15	4
CM703-40	3F x 40 cm Embolectomy Catheter	3	40	0.2	6
CM703-80	3F x 80 cm Embolectomy Catheter	3	80	0.2	6
CM704-40	4F x 40 cm Embolectomy Catheter	4	40	0.7	8
CM704-80	4F x 80 cm Embolectomy Catheter	4	80	0.7	8
CM705-40	5F x 40 cm Embolectomy Catheter	5	40	1.5	10
CM705-80	5F x 80 cm Embolectomy Catheter	5	80	1.5	10
CM706-40	6F x 40 cm Embolectomy Catheter	6	40	2.0	12
CM706-80	6F x 80 cm Embolectomy Catheter	6	80	2.0	12
CM707-40	7F x 40 cm Embolectomy Catheter	7	40	2.5	14
CM707-80	7F x 80 cm Embolectomy Catheter	7	80	2.5	14

balloon contains latex

STERILE R

Quality Management System

Registered to ISO 13485:2003 under CMDCAS
Annex V of MDD 93/42/EEC, FDA

PRODUCT LISTING, DESCRIPTION, AND/OR SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE



CardioMed Supplies Inc.

199 Saint David Street
Lindsay, Ontario, Canada
K9V 5K7
www.cardiomed.com

Tel: 705.328.2518
Fax: 705.328.9747
Toll-free: 800.387.9757
mail@cardiomed.com

CC_EmbolectomyCath_international 2013/09/24