

# Perfusion packs

## Circuits made to order

CardioMed heart-lung packs embody the best in quality workmanship and materials. Designed to provide perfusionists with dependable components for open-heart surgery, our packs offer a fully contained, packaged solution that can be customized to meet clinical requirements.

Each heart-lung pack is packaged in a convenient double-wrapped tray with individually wrapped arterial/venous circuit tubing. Placed in a tray and double or single-wrapped (according to your needs), each extracorporeal circuit offers a choice of tubing also separately wrapped and identified.

CardioMed heart-lung packs and accessory lines are fabricated with the highest grade PVC tubing. If preferred, a Tygon tubing option is also available on demand.



**CardioMed**  
Life-sustaining solutions

# Perfusion packs

## ORDERING INFORMATION

CODE	DESCRIPTION	SPECIFICATIONS	PACKAGING
01-24-92	Cardioplegia crystalloid delivery system	With stainless steel coil, biowrapped, with patient delivery lines 4:1 ratio	Twelve (12) systems per box
01-22-95B	Bridgeable Cardioplegia crystalloid delivery system	With PVC coil, biowrapped, with patient delivery lines 4:1 ratio	Twelve (12) systems per box
01-25-89L	Cardioplegia crystalloid delivery system	With PVC coil, biowrapped, complete in plastic bucket with patient delivery lines Available in 1:1, 2:1, 4:1 or 9:1 ratios	Eighteen (18) systems per box
12-20-93LP	Pediatric Cardioplegia crystalloid delivery system	With PVC coil, biowrapped, complete in plastic bucket with patient delivery lines	Eighteen (18) systems per box
17-051	Quick prime line	Simple	Twenty five (25) systems per box
17-054	Quick prime line	Bifurcated	Twenty five (25) systems per box
17-760	Unidirectional purge line	Purge	Fifty (50) systems per box
17-75x	Pressure monitor separator line	For various monitors	Twenty five (25) systems per box

Except where otherwise noted, all products are latex-free

Except where otherwise indicated, products listed are ETO sterilized

Consult your CardioMed representative for more information about cooling jacket and perfusion accessories.

Pack components include but are not limited to table pack lines, arterial and venous circuits, perfusion lines or arterial filter loops.



CardioMed Supplies Inc.  
 199 Saint David Street,  
 Lindsay, ON K9V 5K7 Canada  
 Tel: 705.328.2518 Toll-free: 800.387.9757  
 Fax: 705.328.9747  
 Web: [www.cardiomed.com](http://www.cardiomed.com) Email: [mail@cardiomed.com](mailto:mail@cardiomed.com)

Quality Management System  
 Registered to ISO 13485: 2003 under CMDCAS  
 Annex V of MDD 93/42/EEC, FDA