





# Ultrafilter (Hemoconcentrator Kits)

A CardioMed cardiovascular product

## Ultrafiltration and flow circuitry

CardioMed's **Ultrafilter (Hemoconcentrator) Kit** combines CardioMed's broad experience in building cardiovascular/perfusion kits with a high-quality ultrafiltration membrane.



- ✓ Custom kits are available
- ✓ Quick connect Hanson connectors on all Ultrafilter and Dialysis lines
- ✓ Rapid ultrafiltration rates virtually eliminate the need for a vacuum source
- ✓ Biocompatible Polysulphone Membrane ensures minimal complement activation with no evident synthesis of cytokines
- ✓ No rinse membrane ensures rapid setup times
- ✓ Kit sizes are available from pediatric up to large adult in Ultrafiltrate and Dialysate configurations

## CardioMed hemoconcentrator kits

### ORDERING INFORMATION

CODE	DESCRIPTION	MEMBRANE
CM9400	Ultrafilter (Hemoconcentrator) tubing set	
CM-9400-50E	Ultrafilter (Hemoconcentrator) kit with CM9400 & Ultrafilter	50 Micron
CM-9400-160	Ultrafilter (Hemoconcentrator) kit with CM9400 & Ultrafilter	60 Micron
CM9405	Ultrafilter (Hemoconcentrator) tubing set	
CM-9405-50E	Ultrafilter (Hemoconcentrator) kit with CM9405 & Ultrafilter	50 Micron
CM-9405-160	Ultrafilter (Hemoconcentrator) kit with CM9405 & Ultrafilter	60 Micron
CM9504	Ultrafilter (Hemoconcentrator) tubing set w/ 4L effluent bag	
CM-9504-50E	Ultrafilter (Hemoconcentrator) kit with CM9504 & Ultrafilter	50 Micron
CM-9504-160	Ultrafilter (Hemoconcentrator) kit with CM9504 & Ultrafilter	60 Micron

**All Ultrafilter kits are ETO sterilized and packaged 20 kits per box**

The 9400 series Ultrafilter kits include :  
an Ultrafilter tubing kit with a female posi-lock

The 9405 series kits include :  
a concentrator tubing kit with luer locks and clamps

The 9504 series kits include :  
a concentrator tubing kit with luer locks and clamps and 4L effluent bag

All Ultrafilter kits come with:  
1 – arterial line  
1 – venous line  
1 – ultrafiltrate line with Hanson connector.

All products listed here are latex-free.

Other accessories and filter options are available.  
Contact your sales representative or Customer Service for more details



### CardioMed Supplies Inc.

199 Saint David Street,  
Lindsay, ON K9V 5K7

Tel: 705.328.2518

Toll-free: 800.387.9757

Fax: 705.328.9747

Web: [www.cardiomed.com](http://www.cardiomed.com)

Email: [mail@cardiomed.com](mailto:mail@cardiomed.com)

Quality Management System  
Registered to ISO 13485: 2003 under CMDCAS  
Annex V of MDD 93/42/EEC, FDA