

Trach Pro closed suction catheters

Combining safety and comfort



Latex Free
DEHP Free

 **CardioMed**

Trach Pro closed suction catheters

ORDERING INFORMATION

ITEM	DESCRIPTION	OUTER DIAMETER	INNER DIAMETER	LENGTH	COLOUR BAND
CM-3106	6Fr Trach Pro CSC full catheter kit	2.0 mm	1.4 mm	300 mm	pea green
CM-3108	8Fr Trach Pro CSC full catheter kit	2.7 mm	1.8 mm	300 mm	blue
CM-3110	10Fr Trach Pro CSC full catheter kit	3.3 mm	2.4 mm	300 mm	black
CM-6110	10Fr Trach Pro CSC full catheter kit	3.3 mm	2.4 mm	600 mm	black
CM-3112	12Fr Trach Pro CSC full catheter kit	4.0 mm	2.8 mm	300 mm	white
CM-6112	12Fr Trach Pro CSC full catheter kit	4.0 mm	2.8 mm	600 mm	white
CM-3114	14Fr Trach Pro CSC full catheter kit	4.7 mm	3.2 mm	300 mm	green
CM-6114	14Fr Trach Pro CSC full catheter kit	4.7 mm	3.2 mm	600 mm	green
CM-3116	16Fr Trach Pro CSC full catheter kit	5.3 mm	3.8 mm	300 mm	orange
CM-6116	16Fr Trach Pro CSC full catheter kit	5.3 mm	3.8 mm	600 mm	orange
CM-3206	6Fr Trach Pro CSC catheter assembly with suction valve	2.0 mm	1.4 mm	300 mm	pea green
CM-3208	8Fr Trach Pro CSC catheter assembly with suction valve	2.7 mm	1.8 mm	300 mm	blue
CM-3210	10Fr Trach Pro CSC catheter assembly with suction valve	3.3 mm	2.4 mm	300 mm	black
CM-6210	10Fr Trach Pro CSC catheter assembly with suction valve	3.3 mm	2.4 mm	600 mm	black
CM-3212	12Fr Trach Pro CSC catheter assembly with suction valve	4.0 mm	2.8 mm	300 mm	white
CM-6212	12Fr Trach Pro CSC catheter assembly with suction valve	4.0 mm	2.8 mm	600 mm	white
CM-3214	14Fr Trach Pro CSC catheter assembly with suction valve	4.7 mm	3.2 mm	300 mm	green
CM-6214	14Fr Trach Pro CSC catheter assembly with suction valve	4.7 mm	3.2 mm	600 mm	green
CM-3216	16Fr Trach Pro CSC catheter assembly with suction valve	5.3 mm	3.8 mm	300 mm	orange
CM-6216	16Fr Trach Pro CSC catheter assembly with suction valve	5.3 mm	3.8 mm	600 mm	orange
CM-6SV	6Fr Trach Pro CSC suction valve connector				

Except where otherwise noted, all products are latex-free

Unless otherwise indicated, products listed are sterilized and packaged 10 per box

The Trach Pro closed suction catheter minimizes lower respiratory tract infections and nosocomial cross-infections, utilizes a blunt tip design to protect the epithelial layer of the tracheal wall, maximizes secretion removal through the use of large, staggered side holes, and features 360 degree swivel adaptors to reduce strain on endotracheal and tracheostomy tubes.

Trach Pro closed suction full catheter kits include:

- T connector with sliding on-off valve
- Closed suction catheter with flush port
- Suction valve connector
- Corrugated Tube

Trach Pro catheter assembly with suction valve includes:

- Closed suction catheter with flush port
- Suction valve connector



CardioMed Supplies Inc.

199 Saint David Street, Lindsay, ON K9V 5K7 Canada

Tel: 705.328.2518 Toll-free: 800.387.9757

Fax: 705.328.9747

Web: www.cardiomed.com Email: mail@cardiomed.com

Quality Management System

Registered to ISO 13485: 2003 under CMDCAS

Annex V of MDD 93/42/EEC, FDA