




HSG / SHG CATHETER

A NEW LEVEL OF CLARITY AND CONVENIENCE IN HSG / SHG PROCEDURES

HSG CATHETERS ARE LATEX-FREE BALLOON CATHETERS DESIGNED TO SEAL THE UTERUS DURING INJECTION AND PROVIDE A CLEAR PICTURE OF THE UTERINE CAVITY AND FALLOPIAN TUBES, WITH MINIMAL OBSTRUCTION OF THE INTERNAL OS.



HSG / SHG Catheter

ORDER CODES

Product Code	Description	Box Quantity
CM-COSF	Cervical OS Finder	25
CM-USM	Uterine Sound w/ 11cm Measurements	25
CM-5NSHG	5F Sonohysterography (SHG) Catheter, Balloonless	25
CM-5BHSG1	5F Hysterosalpingography (HSG) Balloon Catheter w/ 1cc Syringe	10
CM-5BHSGT1	Hysterosalpingography (HSG) Tray, 5F Catheter w/ 1cc Syringe	12
CM-HSGK	Hysterosalpingography (HSG) Kit	12

KIT CONTENTS

Welch-Allyn Speculum - Medium
 Lubricating Jelly
 3x Povidone Iodine Ointment
 2x 6" Applicator
 20CC Syringe
 Sanitary Napkin
 18G x 1.5" Needle
 Gauze Pads
 1x 6" Sponge Tipped Swab Stick
 Uterine Sound
 Cervical OS Finder
 2x 30ml Povidone Iodine Solution
 18" x 26" Underpad
 20" Extension Tubing Set
 Inner Tray

TRAY CONTENTS

5F Balloon Catheter w/ 1cc Syringe
 Welch-Allyn Speculum - Medium
 Lubricating Jelly
 2x Povidone Iodine Ointment
 3x 6" Applicator
 20CC Syringe
 Sanitary Napkin
 Gauze Pads
 18" x 26" Underpad
 Inner Tray



199 St. David Street
 Lindsay, Ontario, Canada
 K9V 5K7

1-800-387-9757

Learn more at:

www.CardioMed.com

